

11 Instruments

11.1 Cutting Instruments

11.2 Utility Instruments

11.3 Band Removers

11.4 Debonding Instruments

11.5 Typing Pliers

11.6 Hand Instruments

11.7 Plier Organizers

11.1 Cutting Instruments

1 Pin and Ligature Cutter OIC-100



2 Hard Wire Cutter Black Oxide OIC-101



3 Mini Pin and Ligature Cutter OIC-102



4 Heavy Duty Cutter OIC-103



5 Distal End Cutter OIC-105



11.2 Utility Instruments

1 Bird Beak Plier OIP-100



2 Bird Beak Plier Black

OIP-101



3 Light Wire Plier

OIP-102



4 Bird Beak Plier with Cutter

OIP-103



5 Straight Utility Plier(How Style)

OIP-104



6 Straight Utility Plier(How Style, Small Tips) OIP-105



7 Utility Plier (Weingart Style, Slim Tips) OIP-106



8 Three-Jaw Plier(Straight) OIP-107



9 Jarabak Plier

OIP-108



10 Kim Plier

OIP-109



11 Loop Forming Plier(Tweed Style)

OIP-110



12 Stop (V-Bend) Plier

OIP-111



13 Nance Loop Forming Plier

OIP-112



14 Hollow Chop

OIP-113



15 Cinch Back Plier


OIP-114



16 Detailing Step Plier OIP-115



17 Distal Ending Bender NiTi OIP-116



18 Ligature Wire Forming Plier OIP-117



19 Hollow Chop Arch Forming Plier OIP-118



20 Lingual Arch Forming Plier OIP-119



21 Optical Style Plier OIP-120



22 Utility Plier with Cutter

OIP-121



11.3 Band Removers

1 Posterior Band Remover

OIR-100



2 Posterior Band Remover Enhanced

OIR-101



3 Convertible Cap Remover

OIP-102



- Safely and easily removes convertible caps.
- Instrument places very little force on the tube, making it particularly effective with direct bond tubes

11.4 Debonding Instruments

1 Direct Bond Bracket Remover

OID-100



2 Direct Bond Bracket Remover Posterior

OID-101



3 Adhesive Remover


OID-102



11.5 Typing Pliers

1 Force Module Separator

OIF-100



2 Surgical Hook Crimping Plier

OIF-101



- Two grooved inserts maintain a firm grip on the tube as it is delivered the arch wire.
- Securely crimps the tube at three points to prevent the ends from flaring.
- Lower beaks are thinner to accommodate curvature of arch hook

3 Sliding Hook Crimping Plier OIF-102



4 Rectangular Arch Forming OIF-103



5 Rectangular Arch Forming Plier with one side OIF-104



6 Arch Forming Gauge

OIF-105



11.6 Hand Instruments

1 Direct Bond Tweezer HID-100



- For use with any bracket from a wide twin to a single
- Reverse action hold bracket firmly and releases smoothly

11.7 Plier Organizers

1 Uprighting Plier Organizer



2 Plier Organizer simple style

