

4. Bands

4.1 Loosed Molar Bands

4.2 Loosed Molar Bands kits

4.3 Prewelded Molar Bands

4.4 Molar Bands Tray

4.1 Loosed Molar Bands

FreeFit Molar Bands

- Permanent **Laser Identification** for easy identification
- **Micro-Etched Interior Surface** ensures maximum retention
- **Anatomically Contoured** made it to be fast placement and accurate fit
- **Ideally Tempered** for easy contouring
- **Rolled Edges** made no sharp edges to maximum patient comfort
- **Straight Interproximal Wall** makes bands seating easy
- **Free Prewelding service** including Buccal Tubes, Buccal Attachments, Lingual Attachments.



Package	No.29--44+#	1st Molar	No.29-44+#	2nd Molar
10/pk	UR	B600-06UR1	UR	B600-07UR1
10/pk	UL	B600-06UL1	UL	B600-07UL1
10/pk	LR	B600-06LR1	LR	B600-07LR1
10/pk	LL	B600-06LL1	LL	B600-07LL1
4/pk	UR/UL/LR/LL	B600-0604	UR/UL/LR/LL	B600-0704

4.2 Loosed Molar Bands kits Kits--100/300/600kit

·Starter Kit – 100 Kit B600-6100

[illegible]

2 nd Molar Band Size		29	29+	30	30+	31	31+	32	32+	33	33+	34	34+	35	35+	36	36+	37	37+	38	38+	39	39+	40	40+	41	41+	42	42+	43	43+	44	44+	Total
U/L	R							1	1	2	2	2	2	2	2	2	2	2	2	1	1	1												25
	L							1	1	2	2	2	2	2	2	2	2	2	2	1	1	1												25

Professional Kit – 300 Kit B600-6300

1 st Molar		29	29+	30	30+	31	31+	32	32+	33	33+	34	34+	35	35+	36	36+	37	37+	38	38+	39	39+	40	40+	41	41+	42	42+	43	43+	44	44+	Total
Band	Size																																	
U/L	R				1	2	2	2	2	3	3	4	4	5	5	5	5	5	5	4	4	3	2	2	1	1	1	1	1	1	1			75
	L				1	2	2	2	2	3	3	4	4	5	5	5	5	5	5	4	4	3	2	2	1	1	1	1	1	1	1			75
2 nd Molar		29	29+	30	30+	31	31+	32	32+	33	33+	34	34+	35	35+	36	36+	37	37+	38	38+	39	39+	40	40+	41	41+	42	42+	43	43+	44	44+	Total
Band	Size																																	
U/L	R				1	2	2	2	2	3	3	4	4	5	5	5	5	5	5	4	4	3	2	2	1	1	1	1	1	1	1			75
	L				1	2	2	2	2	3	3	4	4	5	5	5	5	5	5	4	4	3	2	2	1	1	1	1	1	1	1			75

Master Kit -- 600 Kit B600-6600

1 st Molar Band Size		29	29+	30	30+	31	31+	32	32+	33	33+	34	34+	35	35+	36	36+	37	37+	38	38+	39	39+	40	40+	41	41+	42	42+	43	43+	44	44+	Total
U/L	R	1	1	2	2	2	3	4	4	5	5	6	8	8	8	8	8	8	8	8	8	8	6	5	5	4	4	3	2	2	2	1	1	150
	L	1	1	2	2	2	3	4	4	5	5	6	8	8	8	8	8	8	8	8	8	8	6	5	5	4	4	3	2	2	2	1	1	150
2 nd Molar Band Size		29	29+	30	30+	31	31+	32	32+	33	33+	34	34+	35	35+	36	36+	37	37+	38	38+	39	39+	40	40+	41	41+	42	42+	43	43+	44	44+	Total
U/L	R	1	1	2	2	2	3	4	4	5	5	6	8	8	8	8	8	8	8	8	8	8	6	5	5	4	4	3	2	2	2	1	1	150
	L	1	1	2	2	2	3	4	4	5	5	6	8	8	8	8	8	8	8	8	8	8	6	5	5	4	4	3	2	2	2	1	1	150

4.3 Prewelded Molar Bands

Prewelded with all size of MIM Free-Fit serious Buccal Tubes, attachments as you need.



NO.	0.018								0.022			
31	System		Toq.	Rot.	Non		Con		Non		Con	
31+	Maxillary				R		L		R		L	
32	Edge:1st molar		0	10	B131-16UR		B131-16UL		B141-16UR		B241-16UL	
32+			0	10	B130-17UR		B130-17UL		N/A		N/A	
					B131-26UR		B131-26UL		B141-26UR		B241-26UL	
					B130-27UR		B130-27UL		N/A		N/A	

		2nd molar										
33	Sgl	Roth:1st molar	-14	10	B231-16UR	B231-16UL	B241-16UR	B241-16UL	B231-26UR	B231-26UL	B241-26UR	B241-26UL
33+		2nd molar	-14	10	B230-17UR	B230-17UL	N/A	N/A	B230-27UR	B230-27UL	N/A	N/A
34		Strt Wire	-20	4	B331-16UR	B331-16UL	N/A	N/A	B331-26UR	B331-26UL	N/A	N/A
34+												
35	Dbl (tube with	Edge	-14	10	B134-16UR	B134-16UL	B144-16UR	B144-16UL	B134-26UR	B134-26UL	B144-26UR	B144-26UL
35+	Occlusal	Roth	-25	4	B234-16UR	B234-16UL	B244-16UR	B244-16UL	B234-26UR	B234-26UL	B244-26UR	B244-26UL
36	headgear .004 5)	Strt Wire	-20	4	B334-16UR	B334-16UL	N/A	N/A	B334-26UR	B334-26UL	N/A	N/A
36+												
37	Dbl (Two	Edge	0	10	B132-16UR	B132-16UL	B142-16UR	B142-16UL	B132-26UR	B132-26UL	B142-26UR	B142-26UL
37+	rectangular tube)	Roth	-25	4	B232-16UR	B232-16UL	B242-16UR	B242-16UL	B232-26UR	B232-26UL	B242-26UR	B242-26UL
38												
38+	Tpl	Edge	-14	10	B143-16UR	B143-16UL	N/A	N/A	B143-26UR	B143-26UL	N/A	N/A
39		Roth	-25	4	B243-16UR	B243-16UL	N/A	N/A	B243-26UR	B243-26UL	N/A	N/A
39+												
40	Mandibulary				R	L	R	L	R	L	R	L
40+		Edge:1st molar	0	4	B131-16LR	B131-16LL	B141-16LR	B241-16LL	B131-26LR	B131-26LL	B141-26LR	B241-26LL
41		2nd molar	0	4	B130-17LR	B130-17LL	N/A	N/A	B130-27LR	B130-27LL	N/A	N/A
41+	Sgl	Roth:1st molar	-25	4	B231-16LR	B231-16LL	B241-16LR	B241-16LL	B231-26LR	B231-26LL	B241-26LR	B241-26LL
42		2nd molar	-25	4	B230-17LR	B230-17LL	N/A	N/A	B230-27LR	B230-27LL	N/A	N/A
42+		Strt Wire	-20	4	B331-16LR	B331-16LL	N/A	N/A	B331-26LR	B331-26LL	N/A	N/A
43												
43+	Dbl (Two	Edge	0	4	B132-16LR	B132-16LL	B142-16LR	B142-16LL	B132-26LR	B132-26LL	B142-26LR	B142-26LL
44	rectangular tube)	Roth	-25	4	B232-16LR	B232-16LL	B242-16LR	B242-16LL	B232-26LR	B232-26LL	B242-26LR	B242-26LL

4.4 Molar Bands Trays B601-0600

Easily store and organize your bands and tubes assemblies!

